

Knowledgebase > Using the Care Tracker > How can I share tracking access with other relatives?

How can I share tracking access with other relatives?

Mari Opstad - 2025-11-26 - Comments (0) - Using the Care Tracker

How can I share tracking access with other relatives?

In short:

You can give other relatives access to follow the watch in the **Care Tracker app**. This makes it easy for multiple family members to stay informed and help ensure the user's safety.

Here's how:

- 1. Open the Care Tracker app.
- 2. Select the watch you want to share access to, and tap Settings.
- 3. Choose Alarmreceivers.
- 4. Tap Add alarmreceiver.
- 5. Enter the **name**, **phone number**, **and email address** of the person you want to grant access to. The person will receive an email with a link to create their account in the Care Tracker app.
- 6. Once the new user has registered, they'll be able to see the watch's location and receive alerts just like you.

Tip:

You can remove or change access at any time under **Alarmreceivers**.

Note:

Only the people you actively add as alert recipients can view the watch's location and receive notifications. No one else has access to this information.